



X⁶ Deluxe Chair Bed

The Care Quip X⁶ Deluxe Chair Bed has been designed for the treatment and prevention of pressure problems related to long-term seating.

Using a contoured dual layer air cushion overlaid with pressure relieving high density foam, the X⁶ Chair Bed evenly distributes the occupant's weight, ensuring that there are no main pressure points that can restrict circulation.

The Care Quip X⁶ Deluxe Chair Bed is a must have in facility prevention and treatment of pressure sores. This cost effective, mobile chair bed improves facility care and equips carers with quality pressure care options.



Seating Surface

Four way Stretch Vapour permeable fabric Seating surface increases patient comfort and is designed to improve patient pressure care.

Vapour permeable - water resistant yet breathable and improves skin microclimate

Belgium Made fabric meets European fire standards for being Halogen free fire retardancy

All covers are removable for infection control management

Push Handle

Full Hand tilt and Recline levers allows seamless Patient comfort

Vertical and Horizontal Push Handle positions increases movability and chair control





Safe Handling

Three individual locking castors, along with one directional castor assists with manoeuvring, and the ergonomic handles assist with leverage

Swing Away Backrest Laterals and Armrests

carers demands when positioning the patient

Allows for easy patient side transfers

Specifications

Seat Height	650mm	Overall Depth	1130mm
Seat Depth	560mm	Recliner Length	1790mm
Seat Width	500mm	Tilt Angle	25°
Seat to Top of Back	730mm	Leg Raise Angle	90°
Overall Height	1290mm	Recliner Angle	50°
Overall Width	700mm	S.W.L	180KG

Product Features



Shaped foam and air chambers promotes exceptional pressure redistribution

Air flow chambers helps reduce skin moisture when seated for extended periods of time

High Density foam outer supports improves back posture and redistributes pressure

Foam and air chambers can be removed and customised to the patient's pressure requirements



Foot Plates

Removable comfort padded foot plates to aid patient comfort and heal pressure remediation

Increases patient access, and decreases

Air Seat Section

Contoured Adjustable air cushion increases patient emersion and pressure care

Air Cushion is overlayed with a high-density pressure relieving foam that increases patient stability and comfort